

KneuraSense: Complete Setup s User Guide

1. How to Wear the Device Correctly

For the KneuraSense to accurately track your movements, the knee support must be worn the right way.

- **Positioning:** Wrap the knee support around your leg. Make sure the sensors are centered and positioned directly over your knee joint.
- **Secure Fit:** Fasten the velcro straps firmly. The sensors need to maintain constant contact with your skin to get accurate heart rate and skin temperature readings.

2. Wi-Fi Setup

Your device needs to connect to your home internet to fetch weather updates and securely save your daily health data.

1. **Power On:** Turn on the main power switch. The light will turn **Solid Red** while the device runs a quick self-check.
2. **Setup Mode:** If the device doesn't recognize a Wi-Fi network, the light will **Solid Yellow** and the device will vibrate, letting you know it is ready to be set up.
3. **Connect:** Open the Wi-Fi settings on your smartphone and connect to the network named **SETUP-KneuraSense-DEVICE**.
4. **Enter Password:** A setup screen will automatically pop up on your phone. Simply select your home Wi-Fi name and enter your password.
5. **Success:** Once connected, the device will vibrate to let you know it was successful. It will then restart and securely connect to the system.

3. Account Creation s Device Registration

Before the device can send your health data, it needs to be linked to your personal profile using its unique Device ID.

1. **Locate the ID:** Find the 12-character code (for example: A1:B2:C3:D4:E5:F6) printed on the back label of your KneuraSense central pod.
2. **Web Setup:** Go to the [kneura-sense-koa.vercel.app](#) and create your account. Once logged in, go to the **Settings** menu and type in your **device's ID** to link it to your account.
3. **Clinician Setup (Alternative):** If your doctor is setting this up for you, just give them this ID code and your email address. They will register the device, and you will get an email inviting you to view your dashboard.

4. Battery s Charging

To make sure the device monitors you without interruption, keep it charged.

- **Charging:** Plug the device into a power source. The light will turn **Solid Red** to show that it is currently charging.
- **Fully Charged:** When the battery is completely full, the light will turn **Solid Blue**, meaning it is ready to use.

5. Standing Baseline Calibration

Every time you put the sleeve on, the device needs 5 seconds to learn your natural standing posture so it can measure your movements accurately.

1. **Posture:** Stand perfectly straight with your weight evenly distributed on both legs.
2. **Start:** Press the physical Calibration Button on the device.
3. **Hold Still:** Do not move at all for 5 seconds. The device is taking quick measurements to find your baseline. If you shift your weight or if the sensors aren't strapped tightly enough to your skin, the process will cancel itself to prevent bad readings, and you will need to try again.

6. Interpreting Real-Time Alerts

The smart chip inside the device constantly analyzes your movements and checks your stress score against your personalized safe limits. It will give you live feedback:

- **Solid Green (Normal):** The system is working, and your current joint stress is safe.
- **Solid Yellow (Warning):** Your stress score is getting high (70-85% of your limit). This is a gentle reminder to slow down or correct your posture.
- **Solid Red + Vibration (Risk):** Your joint load is too high for your current incline context and the weather. The device will vibrate continuously. Stop what you are doing and sit down to rest.
- **Flashing Red (Error):** A hardware or system error happened. Turn the device off and back on again.

7. Offline Failsafe Operation

If you go for a walk outside your Wi-Fi range or lose your mobile hotspot connection, don't worry. The device will automatically switch to "Offline Mode" and keep monitoring you safely. **Please connect back to a Wi-Fi network within 1 hour** so the device can successfully save all of your recent activity to your dashboard.